

# Murphy Ring Splint® from DanMic Global

To buy the Murphy ring splint on-line, please follow this link:

[www.danmicglobal.com/murphyringsplint.aspx](http://www.danmicglobal.com/murphyringsplint.aspx)



Murphy Ring Splint® Orthosis

## Sizing

Your Medical Professional can help you determine the proper size using the Murphy ring sizer or a standard ring sizer to measure the involved joint. For half sizes, use the next larger size. Alternative sizing method is to measure the circumference of the joint in centimeters (reference chart below)

Ring Size	Circumference of Joint
2	4.1cm
3	4.4cm
4	4.7cm
5	5.0cm
6	5.3cm
7	5.7cm
8	6.0cm
9	6.3cm
10	6.6cm
11	6.9cm
12	7.2cm
13	7.5cm



DanMic Global - Serving your worldwide therapy product needs  
2059 Camden Avenue #281  
San Jose • CA • 95124